

Safety round-shape edge for subgingival




"Thin and flexible edge part like Explorer"
 Soft insertion with much less damage to the periodontal tissues.

"Flat-surface inside the edge"
 This flat shape fit on the tooth surface and catch the biofilm.

"Round-shape cutting edge with delicacy"
 Prevent from over-treatment and damaging prosthesis

"Same angle with Gracey Curette #11-12"
 Easy control at subgingival, molar and root part.



Gracey Curette DBI

Pursue the MI, Minimum Internation, with delicacy

DBI = Debridement Instrument

Blade Breaker Cvd

- For ophthalmology (eye surgery) scalpel.
- Simply hold and pick the scalpel. Handle part is with lock function.

Focusing Angled tip part for clear view through the operation.



\$2300/pc



No.32362


Size / 120mm

[How to use]



① Hold the scalpel and confirm handle part is locked. Twist to downside and pick the scalpel.

② Now it's ready to use.



With spring

This spring open the tip and make us easy to control by one hand.